


TASK ORDER #

## Important Instructions for OCH Use



**OCH**  
MW 704

Prepared by Dr. Paul B Savage's lab according to the procedures described in "Effects of Lipid Chain Lengths in alpha-Galactosylceramides on Cytokine Release by Natural Killer T Cells" by Randal D. Goff, Ying Gao, Jochen Mattner, Depeng Zhou, Ning Yin, Carlos Cantu, III, Luc Teyton, Albert Bendelac, and Paul B. Savage, *J Am Chem Soc*, **2004**, *126*, 13602-13603.

This solution should be stored at 4°C.

**Each vial contains 1 mL of a 1 mg / mL DMSO stock solution.**

Reagents are to be used for non-commercial research purposes only and may not be used in humans or for clinical diagnosis.